



Rabbit Anti-Human NEK7 monoclonal antibody, clone I.702.5 (CABT-L1492)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	NEK7
Immunogen	Synthetic peptide corresponding to the carboxy terminus of human NEK7
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	I.702.5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Introduction	NIMA-related kinases share high amino acid sequence identity with the gene product of the <i>Aspergillus nidulans</i> 'never in mitosis A' gene, which controls initiation of mitosis.
--------------	--

Keywords

NEK7;NIMA (never in mitosis gene a)-related kinase 7;serine/threonine-protein kinase
Nek7;nimA-related protein kinase 7;never in mitosis A-related kinase 7

GENE INFORMATION

Entrez Gene ID

[140609](#)

UniProt ID

[Q8TDX7](#)
