



Rabbit Anti-Human Naked2 monoclonal antibody, clone T.003.8 (CABT-L1497)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	NKD2
Immunogen	Synthetic peptide corresponding to residues surrounding Ala217 of human Naked2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	T.003.8
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Introduction	In the mouse, Nkd is a Dishevelled-binding protein that functions as a negative regulator of the Wnt-beta-catenin-Tcf signaling pathway.
---------------------	--

Keywords

NKD1;naked cuticle homolog 1 (Drosophila);Naked1;protein naked cuticle homolog 1;naked-1;naked cuticle-1;Dvl-binding protein

GENE INFORMATION

Entrez Gene ID[85409](#)**UniProt ID**[Q969F2](#)
