



# Rabbit Anti-Human Puromycin N-acetyltransferase monoclonal antibody, clone 2I23M9 (CABT-L1510)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	PAC
<b>Immunogen</b>	Protein corresponding to Human Pac (aa 1-199)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Clone</b>	2I23M9
<b>Purification</b>	Protein A Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Concentration</b>	0.5 mg/ml
<b>Buffer</b>	Bacteria
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

## BACKGROUND

**Introduction**

Puromycin N-acetyl transferase functions to detox the cell of puromycin and plays a role in antibiotic resistance.

---

**Keywords**

PAC

---