



Rabbit Anti-Mouse PASK monoclonal antibody, clone U.782.8 (CABT-L1512)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	PASK
Immunogen	Synthetic peptide corresponding to the sequence of mouse PASK
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	U.782.8
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords	PASK;PAS domain containing serine/threonine kinase;PAS domain-containing serine/threonine-protein kinase;KIAA0135;PASKIN;STK37;per-arnt-sim (PAS) domain
----------	--

GENE INFORMATION

Entrez Gene ID [269224](#)

UniProt ID [Q8CEE6](#)
