



Rabbit Anti-Human PLXNA4 monoclonal antibody, clone G.677.3 (CABT-L1534)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	PLXNA4
Immunogen	GST-fusion protein corresponding to residues from the human sequence of Plexin A4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	G.677.3
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords	PLXNA4;plexin A4;PLEXA4;PLXNA4A;PLXNA4B;FAYV2820;PRO34003;plexin-A4
----------	---

GENE INFORMATION

Entrez Gene ID [91584](#)

UniProt ID [Q9HCM2](#)
