



Rabbit Anti-Human AMPK beta-1,2 monoclonal antibody, clone F.538.7 (CABT-L1544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is not cross-reactive with other related proteins.
Target	PRKAB2
Immunogen	Synthetic peptide corresponding to the sequence surrounding His233 residues of human AMPK β 1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Hamster, Human, Mouse, Non-human primate, Rat
Clone	F.538.7
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	FC, ICC, IHC-P, IF, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100 μ g/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords

PRKAB2;protein kinase, AMP-activated, beta 2 non-catalytic subunit;5-AMP-activated protein kinase subunit beta-2;AMPK beta 2;AMPK beta-2 chain;AMPK subunit beta-2;5-AMP-activated protein kinase, beta-2 subunit;AMP-activated protein kinase beta 2 non-catalytic subunit

GENE INFORMATION

Entrez Gene ID[5565](#)**UniProt ID**[O43741](#)
