



# Rabbit Anti-Human Radixin monoclonal antibody, clone L.464.0 (CABT-L1564)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody is not cross-reactive with ezrin or moesin.
<b>Target</b>	RDX
<b>Immunogen</b>	Synthetic peptide corresponding to residues near the carboxy terminus of human radixin
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Non-human primate, Rat
<b>Clone</b>	L.464.0
<b>Purification</b>	Affinity Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Buffer</b>	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
<b>Preservative</b>	See individual product datasheet
<b>Storage</b>	-20°C

## BACKGROUND

**Keywords**RDX;radixin;DFNB24

---

## GENE INFORMATION

**Entrez Gene ID**[5962](#)

---

**UniProt ID**[P35241](#)

---