



Rabbit Anti-Human RSK1/RSK2/RSK3 monoclonal antibody, clone G.050.8 (CABT- L1576)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	RPS6KA1; RPS6KA3; RPS6KA2
Immunogen	Synthetic peptide corresponding to residues surrounding Asp659 of human p90RSK
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Guinea pig, Human, Mouse, Non-human primate, Rat
Clone	G.050.8
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ICC, IP, IF, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords

CLS;D17Wsu134e;HU-1;HU-2;HU-3;ISPK-1;ISPK1;Mapkapk-1a;MAPKAPK-1b;MAPKAPK1A;MAPKAPK1B;MAPKAPK1C;MPK-9;MRX19;p90-Rsk1;p90-RSK2;p90-RSK3;p90Rsk;p90Rsk1;p90RSK2;p90RSK3;p90S6K;pp90rsk;pp90RSK2;pp90RSK3;RGD1563860;Rps6krs1;RPS6KA1;RPS6KA2;RPS6KA3;RSK;Rsk-1;RSK1;RSK2;RSK3;S6K-alpha;S6K-alpha-1;S6K-alpha2;S6K-alpha3
