



Rabbit Anti-Mouse SCD1 monoclonal antibody, clone H.656.20 (CABT-L1588)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	SCD1
Immunogen	Synthetic peptide corresponding to residues surrounding Leu35 of mouse SCD1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	H.656.20
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ICC, IHC-P, IP, IF, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Introduction	Stearoyl-CoA desaturase.
--------------	--------------------------

Keywords

SCD1;F27J15.17;F27J15_17;SCD1;STOMATAL CYTOKINESIS-DEFECTIVE 1;protein
STOMATAL CYTOKINESIS-DEFECTIVE 1

GENE INFORMATION

Entrez Gene ID[20249](#)**UniProt ID**[P13516](#)