



# Rabbit Anti-STAT4 monoclonal antibody, clone I.67.0 (CABT-L1608)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	STAT4
Immunogen	Synthetic peptide corresponding to residues surrounding Lys151 of Stat4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	I.67.0
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ChIP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

## BACKGROUND

Keywords	STAT4;signal transducer and activator of transcription 4
----------	----------------------------------------------------------