



Mouse Anti-CHIKV E2 monoclonal antibody, clone 16A12 (CABT-L1612M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Mouse Anti CHIKV E2 antibody is specific for the E2 envelope protein of Chikungunya virus. It demonstrates negligible cross-reactivity with other members of the alphavirus family, including Western Eastern and Venezuelan Encephalitis viruses. There is no cross-reactivity to Zika, Dengue or other flavivirus antigens.
Target	Chikungunya virus E2 protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	CHIKV
Clone	16A12
Conjugate	Unconjugated
Applications	ELISA, WB, IF
Format	Liquid
Concentration	Lot specific
Size	200 µg
Buffer	PBS, pH7.4
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice