



Rabbit Anti-Human AS160 monoclonal antibody, clone H.964.8 (CABT-L1625)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	TBC1D4
Immunogen	Synthetic peptide corresponding to residues surrounding Ala195 of human AS160
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	H.964.8
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords	TBC1D4;TBC1 domain family, member 4;AS160;TBC1 domain family member 4;akt substrate of 160 kDa;Acrg embryonic lethality minimal region ortholog;TBC (Tre-2, BUB2, CDC16)
----------	--

GENE INFORMATION

Entrez Gene ID [9882](#)

UniProt ID [O60343](#)
