



Rabbit Anti-Human TCF7 monoclonal antibody, clone H.563.0 (CABT-L1630)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is not cross-reactive with LEF1.
Target	TCF7
Immunogen	Synthetic peptide corresponding to a region surrounding Leu158 of human TCF1 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	H.563.0
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords

TCF7;transcription factor 7 (T-cell specific, HMG-box);transcription factor 7;TCF 1;T-cell-specific transcription factor 1;TCF-1;FLJ36364;MGC47735

GENE INFORMATION

Entrez Gene ID[6932](#)**UniProt ID**[P36402](#)
