



Rabbit Anti-Human TNK1 monoclonal antibody, clone F.460.20 (CABT-L1638)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is not cross-reactive with other related proteins.
Target	TNK1
Immunogen	Synthetic peptide corresponding to the sequence around Pro97 of human TNK1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	F.460.20
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Introduction	TNK1 is a nonreceptor tyrosine kinase. These kinases, like members of the SRC and JAK families, mediate intracellular signaling downstream of receptor activation.
Keywords	TNK1;Tyrosine kinase non receptor 1;non-receptor tyrosine-protein kinase TNK1;CD38 negative kinase 1

GENE INFORMATION

Entrez Gene ID	8711
UniProt ID	Q13470