



Rabbit Anti-Human Tensin 4 monoclonal antibody, clone TQ94 (CABT-L1640)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	TNS4
Immunogen	Synthetic peptide corresponding to C-terminus of human CTEN protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	TQ94
Purification	Protein A/G Purified
Conjugate	Unconjugated
Applications	IHC-P
Positive Control	A suggested positive control is liver hepatocellular carcinoma.
Format	Liquid
Buffer	PBS, pH 7.6, with 1% BSA
Preservative	0.1% Sodium Azide
Storage	4°C, do not freeze

BACKGROUND

Introduction

C-TEN/TNS4 links the focal adhesion and cell adhesion signaling to the actin cytoskeleton. It may be involved in cell survival and migration. The protein may be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton and promote apoptosis, via its cleavage by caspase-3.

GENE INFORMATION

Entrez Gene ID	84951
----------------	-----------------------

UniProt ID	Q8IZW8
------------	------------------------
