



# RHA™ Mouse anti-Chloromycetin monoclonal antibody, clone CAP (CABT-L1668)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Chloromycetin [KLH]
Isotype	IgG
Source/Host	Mouse
Clone	CAP
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C, avoid repeat freeze- thaw cycles