



Rabbit Anti-Human FOXP4 polyclonal antibody (CABT-L1757)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody may react with (Predicted by homology) : Bovine, Dog, Rabbit, Rat
Target	FOXP4
Immunogen	Synthetic peptide corresponding to C-terminus of human FOXP4 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunoaffinity purified
Conjugate	Unconjugated
Applications	IHC-P
Molecular Weight	74 kDa
Cellular Localization	Nucleus
Positive Control	Lung
Format	Liquid
Buffer	PBS, 1% BSA, pH 7.6
Preservative	< 0.1% Sodium Azide
Storage	2-8°C. Do not freeze. The user must validate any other storage conditions. When properly

stored, the reagent is stable to the date indicated on the label. Do not use the reagent beyond the expiration date.

BACKGROUND

Introduction

This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Many members of the forkhead box gene family, including members of subfamily P, have roles in mammalian oncogenesis. This gene may play a role in the development of tumors of the kidney and larynx. Alternative splicing of this gene produces multiple transcript variants, some encoding different isoforms.

Keywords

FOXP4;forkhead box P4;hFKHLA;forkhead box protein P4;winged-helix repressor FOXP4;fork head-related protein like A;fork head-related protein-like A

GENE INFORMATION

Entrez Gene ID

[116113](#)

UniProt ID

[Q9DBY0](#)
