



Rabbit Anti-HSV-1 ICP34.5 (N-terminal) Polyclonal Antibody (CABT-L184N)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	N-terminal residues of HSV-1 ICP34.5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	Lot specific
Size	100 µg
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Herpes simplex virus 1 and 2 (HSV-1 and HSV-2), also known by their taxonomical names Human alphaherpesvirus 1 and Human alphaherpesvirus 2, are two members of the human
--------------	---

Herpesviridae family, a set of new viruses that produce viral infections in the majority of humans. Both HSV-1 (which produces most cold sores) and HSV-2 (which produces most genital herpes) are common and contagious. They can be spread when an infected person begins shedding the virus.

Keywords

HSV;Herpes Simplex Virus;Herpes Simplex Virus 1;HSV-1;ICP34.5;HSV-1 ICP34.5
