



# Mouse Anti-Drosophila Tumbleweed monoclonal antibody, clone 2I6 (CABT-L1938)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Drosophila Tumbleweed
Immunogen	Tumbleweed, full length
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Drosophila melanogaster
Clone	2I6
Conjugate	Unconjugated
Applications	IP, IB
Molecular Weight	69.8 kDa
Format	Liquid, Supernant
Concentration	Lot specific
Size	1 ml
Preservative	None
Storage	Although many cell products are maintained at 4°C for years without loss of activity, shelf-life at 4°C is highly variable. To ensure retention of antibody activity, we recommend aliquotting the product into two parts: 1) a volume of antibody stored at 4°C to be used within two weeks. 2) the remaining product diluted with an equal volume of molecular grade glycerol and stored at -20°C.

# BACKGROUND

Keywords	tum
----------	-----

---

# GENE INFORMATION

Gene Name	tum
-----------	-----

---

Synonyms	tum
----------	-----

---