



Mouse Anti-2,4-D monoclonal antibody, clone 24D1 (CABT-L1951)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated to 2,4-D
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	24D1
Purification	Protein G purified.
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	0.01M TBS(pH7.4) with 1% BSA and 50% Glycerol.
Preservative	None
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords 2,4-Dichlorophenoxyacetic acid;2,4-D

GENE INFORMATION

Gene Name 2,4-Dichlorophenoxyacetic acid

Synonyms 2,4-Dichlorophenoxyacetic acid;2,4-D
