



Mouse Anti-Propoxyphene monoclonal antibody, clone C2789N (CABT-L1966)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	MAb to Propoxyphene. Monoclonal Antibody to Propoxyphene
Specificity	Propoxyphene
Target	Propoxyphene
Immunogen	Propoxyphene [BSA].
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	C2789N
Purification	> 95% pure (SDS-PAGE). Protein A Chromatography
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.02% Sodium Azide

Storage	Short term (up to 7 days) store at 2-8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles
---------	--

BACKGROUND

Keywords	(+)-PROPOXYPHENE; PROPOXYPHENE; (+)-1,2-diphenyl-2-propionoxy-3-methyl-4-dimethylaminobutane; (+)-4-(Dimethylamino)-3-methyl-1,2-diphenyl-2-butanol propionate; (+)-4-dimethylamino-1,2-diphenyl-3-methyl-2-propionyloxybutane; (+)-propoxyphen; (2S,3R)-(+)-4-(Dimethylamino)-3-methyl-1,2-diphenyl-2-butanol propionate(ester); (s-(r*,s*))-(s)-alpha-(2-(dimethylamino)-1-methylethyl)-alpha-phenylbenzeneethanolpropanoa; 1,2-diphenyl-2-propionoxy-3-methyl-4-dimethylaminobutane
----------	--

GENE INFORMATION

Gene Name	(+)-PROPOXYPHENE
Synonyms	(+)-PROPOXYPHENE; PROPOXYPHENE; (+)-1,2-diphenyl-2-propionoxy-3-methyl-4-dimethylaminobutane; (+)-4-(Dimethylamino)-3-methyl-1,2-diphenyl-2-butanol propionate; (+)-4-dimethylamino-1,2-diphenyl-3-methyl-2-propionyloxybutane; (+)-propoxyphen; (2S,3R)-(+)-4-(Dimethylamino)-3-methyl-1,2-diphenyl-2-butanol propionate(ester); (s-(r*,s*))-(s)-alpha-(2-(dimethylamino)-1-methylethyl)-alpha-phenylbenzeneethanolpropanoa; 1,2-diphenyl-2-propionoxy-3-methyl-4-dimethylaminobutane