



# Mouse Anti-Human CD77 monoclonal antibody, clone 6C6 (CABT-L1967)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Human CD77
Immunogen	Human peripheral blood mononuclear cells
Isotype	IgM
Source/Host	Mouse
Species Reactivity	Human
Clone	6C6
Purification	Affinity purification
Conjugate	Unconjugated
Applications	FC, WB, IHC
Format	Liquid
Concentration	Lot specific
Size	100 µg
Buffer	Phosphate-buffered solution, pH 7.2
Preservative	0.09% Sodium Azide
Storage	The antibody solution should be stored undiluted between 2°C and 8°C.

## BACKGROUND

Keywords CD77

---

GENE INFORMATION

Gene Name CD77

---

Synonyms CD77

---