



# Rabbit Anti-Human CYFRA21-1 Polyclonal Antibody (CABT-L2128)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Polyclonal Antibody to Cytokeratin Fragment Antigen 21-1 (Knockout Validated)
<b>Specificity</b>	The antibody is a rabbit polyclonal antibody raised against CYFRA21-1. It has been selected for its ability to recognize CYFRA21-1 in immunohistochemical staining and western blotting.
<b>Target</b>	CYFRA21-1
<b>Immunogen</b>	Recombinant fragment corresponding to human CYFRA21-1 (Asp244~Leu400)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	200 µg
<b>Buffer</b>	Supplied as solution form in 0.01M PBS with 50% glycerol, pH7.4.
<b>Preservative</b>	0.05% Proclin-300

<b>Storage</b>	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
<b>Ship</b>	4°C with ice bags

## BACKGROUND

<b>Introduction</b>	CYFRA21-1 is a commonly used tumor marker in lung cancer diagnosis and research.
<b>Keywords</b>	CYFRA21-1

## GENE INFORMATION

<b>Official Symbol</b>	CYFRA21-1
------------------------	-----------