



Human Anti-AAV9 (intact particle) monoclonal antibody, clone BEL0 (CABT-L21Y197)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Anti-AAV9 (intact particle) monoclonal antibody, clone BEL0 This chimeric IgG1 consists of the human constant region of the heavy chain, mouse variable region of the heavy chain and mouse light chain. The human constant region of the heavy chain can be directly recognized by the anti-human conjugate.
Immunogen	recombinant produced adeno-associated virus type 9 (AAV9), carrying GFP
Isotype	IgG1
Source/Host	Human
Species Reactivity	AAV9
Clone	BEL0
Purification	Affinity purified
Conjugate	Unconjugated
Applications	ELISA, Control Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Size	100 µg, 1 mg
Buffer	PBS with Preservative
Storage	Store at -20°C.

Ship

Wet ice

GENE INFORMATION

Synonyms

AAV9; AAV 9; Adeno-associated virus; Adeno-associated virus 9; Adeno-associated virus type 9
