



Rabbit Anti-Human GPR4 (aa 283-362) polyclonal antibody (CABT-L225R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within human GPR4 aa 283-362. Sequence: PILYCLVNEGARSDVAKALHNLLRFLASDKPQEMANASLTLETPLTSKRNSTAKAMTGS WAATPPSQGDQVQLKMLPPAQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in WB. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Preservative	0.02% sodium azide
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Keywords

GPR4; G protein-coupled receptor 4; G-protein coupled receptor 4; G protein coupled receptor 19; G protein coupled receptor 4; GPR 4; GPR19; GPR4; Probable G protein coupled receptor GPR4; G-protein coupled receptor 19;
