



Anti-Cytomegalovirus chimeric monoclonal antibody, clone A4B4C9 (CABT-L2395)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	It is a Mouse/Human chimeric monoclonal antibody produced in transgenic mice by replacing the mouse sequence of the heavy chain constant region (IgM, IgG or IgA loci) by the corresponding human sequence. After immunization with the antigen of interest, generated antibody clones are cultivated by standard hybridoma techniques. They consist of the human constant region of the heavy chain, mouse variable region of the heavy chain and mouse light chain. The human constant region of the heavy chain can be directly recognized by the anti-human conjugate, which is used in numerous in vitro diagnostic assays.
Specificity	This antibody recognizes Cytomegalovirus (CMV)
Target	Cytomegalovirus
Isotype	IgM
Source/Host	Mouse
Species Reactivity	Virus
Clone	A4B4C9
Purification	Unpurified
Conjugate	Unconjugated
Applications	ELISA
Preparation	The antibody has been generated in transgenic mice whose sequence for the IgM heavy chain constant region is replaced by the corresponding human sequence. After immunization of mice, a hybridoma cell line has been established. The antibody is produced industrially by standard hybridoma cell line techniques under sterile conditions. The antibody is presented in cell culture

supernatant.

Format	Liquid
Size	1 ml
Buffer	This cell culture supernatant is supplied in Iscove's Modified Dulbecco's Medium (IMDM), supplemented with 5% FBS, 1% L-Glutamine, 1% Penicillin/Streptomycin, 50 µM 2-Mercaptoethanol.
Preservative	0.01% Sodium Azide
Storage	2–8 °C. Do not use if turbid.
Ship	Wet ice

BACKGROUND

Introduction Cytomegalovirus is a member of the herpes virus group, which includes herpes simplex virus types 1 and 2, varicella zoster virus (which causes chicken pox), and Epstein Barr virus (which causes infectious mononucleosis). These viruses share a characteristic ability to remain dormant within the body over a long period.

Keywords CMV;Cytomegalovirus;HCMV

GENE INFORMATION

Synonyms CMV; Cytomegalovirus; HCMV