



# Anti-Measles chimeric monoclonal antibody, clone Q6I7 (CABT-L2396)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	It is a Mouse/Human chimeric monoclonal antibody produced in transgenic mice by replacing the mouse sequence of the heavy chain constant region (IgM, IgG or IgA loci) by the corresponding human sequence. After immunization with the antigen of interest, generated antibody clones are cultivated by standard hybridoma techniques. They consist of the human constant region of the heavy chain, mouse variable region of the heavy chain and mouse light chain. The human constant region of the heavy chain can be directly recognized by the anti-human conjugate, which is used in numerous in vitro diagnostic assays.
<b>Specificity</b>	This antibody recognizes Measles
<b>Target</b>	Measles
<b>Isotype</b>	IgM
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Virus
<b>Clone</b>	Q6I7
<b>Purification</b>	Unpurified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Preparation</b>	The antibody has been generated in transgenic mice whose sequence for the IgM heavy chain constant region is replaced by the corresponding human sequence. After immunization of mice, a hybridoma cell line has been established. The antibody is produced industrially by standard hybridoma cell line techniques under sterile conditions. The antibody is presented in cell culture

supernatant.

<b>Format</b>	Liquid
<b>Size</b>	1 ml
<b>Buffer</b>	This cell culture supernatant is supplied in Iscove's Modified Dulbecco's Medium (IMDM), supplemented with 5% FBS, 1% L-Glutamine, 1% Penicillin/Streptomycin, 50 µM 2-Mercaptoethanol.
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	2–8 °C. Do not use if turbid.
<b>Ship</b>	Wet ice

## BACKGROUND

<b>Introduction</b>	Measles (also sometimes known as English measles), also known as morbilli, is an infection of the respiratory system caused by a virus, specifically a paramyxovirus of the genus Morbillivirus.
<b>Keywords</b>	morbilli;

## GENE INFORMATION

<b>Synonyms</b>	morbilli;
-----------------	-----------