



Anti-Rubella chimeric monoclonal antibody, clone B8B22I4 (CABT-L2400)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	It is a Mouse/Human chimeric monoclonal antibody produced in transgenic mice by replacing the mouse sequence of the heavy chain constant region (IgM, IgG or IgA loci) by the corresponding human sequence. After immunization with the antigen of interest, generated antibody clones are cultivated by standard hybridoma techniques. They consist of the human constant region of the heavy chain, mouse variable region of the heavy chain and mouse light chain. The human constant region of the heavy chain can be directly recognized by the anti-human conjugate, which is used in numerous in vitro diagnostic assays.
Specificity	This antibody recognizes Rubella
Target	Rubella
Isotype	IgM
Source/Host	Humanized
Species Reactivity	Virus
Clone	B8B22I4
Purification	Unpurified
Conjugate	Unconjugated
Applications	ELISA
Preparation	The antibody has been generated in transgenic mice whose sequence for the IgM heavy chain constant region is replaced by the corresponding human sequence. After immunization of mice, a hybridoma cell line has been established. The antibody is produced industrially by standard

hybridoma cell line techniques under sterile conditions. The antibody is presented in cell culture supernatant.

Format	Liquid
---------------	--------

Buffer	This cell culture supernatant is supplied in Iscove's Modified Dulbecco's Medium (IMDM), supplemented with 5% FBS, 1% L-Glutamine, 1% Penicillin/Streptomycin, 50 µM 2-Mercaptoethanol.
---------------	---

Preservative	0.09% Sodium Azide
---------------------	--------------------

Storage	2–8 °C. Do not use if turbid.
----------------	-------------------------------

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Rubella virus is classified as a togavirus, genus Rubivirus. It is most closely related to group A arboviruses, such as Eastern and Western Equine Encephalitis viruses. It is an enveloped RNA virus, with a single antigenic type that does not cross react with other members of the togavirus group. Rubella virus is relatively unstable and is inactivated by lipid solvents, trypsin, formalin, ultraviolet light, low pH, heat, and amantadine.
---------------------	---

Keywords	Rubellavirus; Rubella virus; Togaviridae; Rubivirus
-----------------	---

GENE INFORMATION

Synonyms	Rubellavirus; Rubella virus; Togaviridae; Rubivirus
-----------------	---
