



Mouse Anti-Phenylethylamine monoclonal antibody, clone PEA (CABT-L2510)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Phenylethylamine
Immunogen	Phenylethylamine with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PEA
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction Phenethylamine (PEA) is an organic compound, natural monoamine alkaloid, and trace amine which acts as a central nervous system stimulant in humans.

Keywords Phenylethylamine;β-phenethylamine;phenethylamine;PEA
