



# Mouse Anti-Bromobuterol monoclonal antibody, clone BL (CABT-L2511)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Cross-Reactivity: Clenbuterol 0.02ppb
Immunogen	Bromobuterol with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	BL
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction**

Among the  $\beta$ -agonist drugs illegally used as repartitioning agents in food producing animals, clenbuterol (CBL) and bromobuterol (BBL) are frequently recovered during routine control by national authorities.

---

**Keywords**

Brombuterol;Bromobuterol;AC1MJ3D4;

---