



# Mouse Anti-Clorprenaline monoclonal antibody, clone CLP (CABT-L2513)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Clorprenaline
Immunogen	Clorprenaline with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CLP
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction**

clorprenaline (uncountable) A beta-adrenergic agonist.

---

**Keywords**

CLORPRENALINE;3811-25-4;Isoprofenamine;Clorprenalina;Isoprophenamine;

---