



Mouse Anti-Alternariol monoclonal antibody, clone AOH (CABT-L2518)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Alternariol
Immunogen	Alternariol with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	AOH
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Alternariol is a toxic metabolite of Alternaria fungi. It is an important contaminant in cereals and fruits. Alternariol exhibits antifungal and phytotoxic activity. It is reported to inhibit cholinesterase enzymes. It is also a mycoestrogen.alternariol is a full androgen agonis

Keywords

Alternariol;3;7;9-Trihydroxy-1-methyl-6H-benzo[c]chromen-6-one
