



Mouse Anti-Tenuazonic acid monoclonal antibody, clone TEA (CABT-L2519)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Tenuazonic acid
Immunogen	Tenuazonic acid with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TEA
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Tenuazonic acid is a mycotoxin produced by *Alternaria* species. It is a powerful eukaryotic protein synthesis inhibitor

Keywords

Tenuazonic acid;(5S)-3-Acetyl-5-[(2S)-butan-2-yl]-4-hydroxy-2,5-dihydro-1H-pyrrol-2-one
