



Mouse Anti-Estrone monoclonal antibody, clone E1 (CABT-L2522)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Estrone
Immunogen	Estrone with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	E1
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Estrone is one of the three naturally occurring estrogens, the others being estradiol and estriol. Estrone is produced primarily from androstenedione originating from the gonads or the adrenal cortex. In premenopausal women, more than 50% of the estrone is secreted by the ovaries. In prepubertal children, men and non supplemented postmenopausal women the major portion of estrone is derived from peripheral tissue conversion of androstenedione. Interconversion of estrone and estradiol also occurs in peripheral tissue. Bioassay data indicate that the estrogenic action is much less than estradiol. Estrone is a primary estrogenic component of several pharmaceutical preparations, including those containing conjugated and esterified estrogens. In premenopausal women estrone levels generally parallel those of estradiol. After menopause estrone levels increase, possibly due to increased conversion of androstenedione to estrone.

Keywords

Estrone
