



Mouse Anti-Apramycin monoclonal antibody, clone APM (CABT-L2527)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Apramycin
Immunogen	Apramycin with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	APM
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Aminoglycosides are a family of bacterial antibiotics that are used in the treatment of specific bacterial infections. They display a concentration dependent killing action and are active against a wide range of aerobic Gram-negative bacilli. Aminoglycosides are molecules that are comprised of an amino group and a sugar group. They operate by inhibiting the bacteria from producing proteins vital to its growth. More specifically, they bind to the bacterial 30S ribosomal subunit where they prevent the translocation of the peptidyl-tRNA from the A-site to the P-site, subsequently giving rise to a misreading of mRNA resulting in the inhibition of protein synthesis. This consequently results in a disruption to the integrity of the bacterial cell membrane. In addition to their use to prevent bacterial infection, aminoglycosides have been used as growth promoters in food producing animals.

Keywords

Apramycin
