



# Mouse Anti-Doramectin monoclonal antibody, clone DOR (CABT-L2531)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Cross-Reactivity: Enamectin 1ppb; Avermectin 1.5ppb; Ivermectin 4ppb
<b>Immunogen</b>	Doramectin with carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	DOR
<b>Purification</b>	Purified from mouse ascites.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	-20°C. Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction**

Doramectin (Dectomax) is a veterinary drug approved by the Food and Drug Administration (FDA) for the treatment of parasites such as gastrointestinal roundworms, lungworms, eyeworms, grubs, sucking lice and mange mites in cattle.

---

**Keywords**

Doramectine;Doramectinum;Doramectin;

---