



Mouse Anti-4,4'- dinitrocarbanilide monoclonal antibody, clone DNC (CABT-L2532)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	4,4'- dinitrocarbanilide
Immunogen	4,4'- dinitrocarbanilide with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DNC
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction	The active component of the antifertility agent nicarbazin
Keywords	1;3-Bis(4-nitrophenyl)urea;587-90-6;4;4'-Dinitrocarbanilide;Urea;N;N'-bis(4-nitrophenyl)-;4;4'-Dinitrodiphenylurea;N;N'-Di(p-nitrophenyl)urea