



Mouse Anti-Flunixin monoclonal antibody, clone FNX (CABT-L2549)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross-Reactivity: 5-hydroxyflunixin 5 ppb
Immunogen	Flunixin with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	FNX
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction	Flunixin is a non-steroidal anti-inflammatory drug (NSAID), analgesic, and antipyretic used in horses, cattle and pigs. It is often formulated as the meglumine salt.
Keywords	Flunixin;