



Mouse Anti-Sildenafil monoclonal antibody, clone SDN (CABT-L2551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Sildenafil
Immunogen	Sildenafil with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	SDN
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Sildenafil, sold as the brand name Viagra among others, is a medication used to treat erectile dysfunction and pulmonary arterial hypertension. Its effectiveness for treating sexual dysfunction in women has not been demonstrated

Keywords

Sildenafil;
