



Mouse Anti-Butafosfan monoclonal antibody, clone BFF (CABT-L2555)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Butafosfan
Immunogen	Butafosfan with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	BFF
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction An organic phosphorus compound used as a drug in animals.

Keywords Butafo;Butapho;BUTAFOSFAN;Butaphosphan;
