



Mouse Anti-Glycochenodeoxycholic acid monoclonal antibody, clone GCDCA (CABT-L2562)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Glycochenodeoxycholic acid
Immunogen	Glycochenodeoxycholic acid with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	GCDCA
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction	Glycochenodeoxycholic acid is a bile salt formed in the liver from chenodeoxycholic acid and glycine, usually found as the sodium salt. It acts as a detergent to solubilize fats for absorption.
Keywords	Glycochenodeoxycholic acid;GCDCA