



Sheep Anti-DHEA-S monoclonal antibody, clone RS21306 (CABT-L259M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Dehydroepiandrosterone Sulphate
Target	Dehydroepiandrosterone Sulphate
Immunogen	DHEAS-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Clone	RS21306
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Buffer	Supplied in 15 mM potassium phosphate, 0.85% NaCl, pH 7.2
Preservative	0.05% Sodium Azide
Storage	Upon receipt, store at 2-8°C
Ship	Wet ice