



# Rabbit Anti-Human QRICH2 (aa 120-220) polyclonal antibody (CABT-L264R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	QRICH2 (glutamine rich 2) is a protein-coding gene.
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids 120-220 of Human QRICH2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Suitable for use in IHC-P. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C.
<b>Ship</b>	Wet ice

# BACKGROUND

## Keywords

QRICH2; glutamine rich 2; glutamine-rich protein 2

---