



Mouse Anti-Bovine Ganglioside GM1 monoclonal antibody, clone HN2-2 (CABT- L271M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This clone recognizes all vertebrate species.
Immunogen	KLH-conjugated bovine brain GM1
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Bovine, Human, Mouse, Rat, Vertebrates
Clone	GM1
Purification	Unpurified
Conjugate	Unconjugated
Applications	ELISA, IF, IHC, WB
Epitope	GM1
Format	Liquid
Size	1 ml
Preservative	See individual product datasheet
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction	GM1 (monosialotetrahexosylganglioside) the "prototype" ganglioside, is a member of the ganglioseries of gangliosides which contain one sialic acid residue.
Keywords	GM1 ganglioside;Monosialoanglioside;Ganglioside GM1

GENE INFORMATION

Official Symbol	Ganglioside GM1
-----------------	-----------------