



Rabbit Anti-Human MGAT5B polyclonal antibody (CABT-L294R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide derived from human MGAT5B: 101-200
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	IHC Recommended working condition FC: 0.25 ug/10 ⁶ cells Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded
Format	Liquid, Purified
Concentration	Lot specific
Size	50 µl, 100 µl, 200 µl
Buffer	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20°C.
Ship	Wet ice