



Mouse Anti-Thiamethoxam monoclonal antibody, clone TMO (CABT-L3071)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Thiamethoxam with carrier protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TMO
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Thiamethoxam is a systemic insecticide in the class of neonicotinoids. It has a broad spectrum of activity against many types of insects.

Keywords

Thiamethoxam;153719-23-4;Diaclofen;Actara;Cruiser;thiamethoxam
