



Rabbit Anti-vc-PAB-MMAE polyclonal antibody (CABT-L3100)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-vc-PAB-MMAE pAb
Target	vc-PAB-MMAE
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Antibody was produced by ascites and then isolated via Protein A/G chromatography
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	1 mg/ml
Size	100 µl, 500 µl
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	For short term storage, store at 4°C up to 6 months from date of opening or thawing. Long time storage is recommended at -20°C. Avoid repeated freeze-thaw cycles.
Ship	Dry ice

BACKGROUND

Introduction	Val-Cit-PAB-MMAE is a drug-linker conjugate for ADC by using the anti-mitotic agent, monomethyl auristatin E (MMAE), linked via the peptide Val-Cit-PAB.
Keywords	VC-PAB-MMAE;VcMMAE;Valine-citrulline-p-aminobenzoyloxycarbonyl-monomethyl auristatin E;Vedotin;monomethyl auristatin E

GENE INFORMATION

Official Symbol	Val-Cit-PAB-MMAE
-----------------	------------------