



Mouse Anti-Paclitaxel monoclonal antibody, clone 464 (CABT-L3222)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Paclitaxel
Immunogen	Paclitaxel
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	464
Purification	Protein A Purified
Conjugate	Unconjugated
Applications	IA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.09% Sodium Azide
Storage	Long time storage is recommended at -20°C.

BACKGROUND

Introduction

Paclitaxel (PTX) is a chemotherapy medication used to treat a number of types of cancer. This includes ovarian cancer, breast cancer, lung cancer, Kaposi sarcoma, cervical cancer, and pancreatic cancer. It is given by injection into a vein. There is also an albumin-bound formulation

Keywords

Paclitaxel;TAXOL;33069-62-4;Taxol A;Abraxane;Paxene;Paxceed;Plaxicel;
