



Human Anti-Human CD314/NKG2D monoclonal Antibody, clone SFB808 [PE] (CABT-L386H-PE)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG1
Source/Host	Human
Species Reactivity	Human
Clone	SFB808
Purification	Affinity purification. Purity:>95% Determined by SDS-PAGE
Conjugate	PE
Applications	Suitable for use in FC Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS, pH 7.0 Contains no stabilizers or preservatives Endotoxin Level: <2EU/mg (<0.002EU/μg) Determined by LAL gel clotting assay
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Keywords

CD314; CD314 antigen; D12S2489E; Killer cell lectin like receptor subfamily K member 1; Killer cell lectin-like receptor subfamily K member 1; KLR; Klrk1; NK cell receptor D; NK lectin-like receptor; NKG2 D activating NK receptor
