



Rat Anti-Mouse IL7 Monoclonal antibody, clone 315737 (CABT-L4077R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Mouse IL-7
Immunogen	E. coli-derived recombinant mouse IL-7
Isotype	IgG2b
Source/Host	Rat
Species Reactivity	Mouse
Clone	315737
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB
Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Format	Lyophilized
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Keywords

Interleukin 17; IL-17
